

Searching for a BNP standard: glycosylated proBNP as a common calibrator enables improved comparability of BNP immunoassays results



Semenov A.¹, Tamm N.¹, Apple F.², Schulz K.², Love S.², Feygina E.³, Katrukha A.^{1,3}

1 - HyTest Ltd., Turku, Finland; 2 - Department of Laboratory Medicine and Pathology, Hennepin County Medical Center and University of Minnesota, Minnesota, USA; 3-School of Biology, Moscow State University, Moscow, Russia

Introduction

B-type natriuretic peptide (BNP) is widely accepted as a useful and cost-effective biomarker of cardiac function. Clinical studies suggest that absolute plasma BNP concentrations measured by various immunoassays differ substantially, complicating interpretation of results and rendering the cut-off concentration method dependent. One of the reasons for the lack of equivalence between existing BNP immunoassays may be the absence of a common calibrator. Presently, there is no agreement on which BNP or peptide standard should be used for calibration of BNP assays, as manufacturers are using different BNP calibrators. Considering this, we suggested that a common calibrator may reduce the degree of between-assay variability of existing commercial BNP immunoassays. Since endogenous proBNP is O-glycosylated, such common BNP calibrator also is expected to be glycosylated.

THE PURPOSE OF OUR STUDY: to compare BNP-related proteins to determine a form that could be used as a common calibrator to improve the comparability of commercial BNP immunoassays.

Materials and Methods

PLASMA SAMPLES: EDTA-plasma samples were obtained from 20 acute and 20 chronic heart failure patients at Hennepin County Medical Center, Minneapolis, MN. Plasma collection was performed in the presence of protease inhibitors to prevent proteolytic degradation of BNP (benzamidine 10 mmol/L and AEBF 5 mmol/L).

BNP IMMUNOASSAYS: BNP concentrations were measured with five commercial BNP assays: Alere Triage, Siemens Centaur XP, Abbott I-STAT, Beckman Access2 (performed at Hennepin County Medical Center, Minneapolis, MN) and ET Healthcare Pylon BNP assay (performed at ET Healthcare) (Fig. 1).




FIGURE 1. Antibodies and standard materials of BNP immunoassays used in the study (according to the information presented at www.ifcc.org).

EXTERNAL CALIBRATORS: synthetic BNP (Bachem), recombinant BNP expressed in *Escherichia coli* (*E. coli*) (Raybiotech), recombinant nonglycosylated proBNP (expressed in *E. coli*; HyTest), His-tagged (N-terminal) recombinant nonglycosylated proBNP (expressed in *E. coli*; Raybiotech), recombinant glycosylated proBNP (expressed in HEK cells; HyTest) and recombinant glycosylated proBNP (expressed in CHO cells, in-house).

CALCULATION OF THE BNP VALUES WITH EXTERNAL CALIBRATORS: the initial BNP concentrations were obtained with the internal standards provided by a manufacturer. External calibrators were used to create calibration curves, which were plotted in logarithmic scale (log-log). The equations obtained with power-law fitting for every calibrator were used to recalculate the BNP concentrations in plasma samples.

The degree of equivalence was analyzed by the linear regression relationship between results of every pair of BNP immunoassays for internal calibrators and 6 external calibrators.

Results

We observed that the different BNP assays were not equal in recognition of BNP, nonglycosylated proBNP and glycosylated proBNP, reflecting the differences in cross-reactivity of commercial BNP immunoassays for different BNP-related forms (Fig. 2).

Absolute BNP concentrations measured in HF plasma samples differed considerably between BNP immunoassays for internal calibrators (up to 3.6-fold differences; range 0.9 to 3.6) and five out of six external calibrators. When glycosylated proBNP (expressed in HEK cells) was utilized as a common external calibrator for all assays, a significant reduction of the between-assay variability was achieved with regression line slopes close to 1.0 for almost every pair of assays (Table 1).

Two forms of glycosylated proBNP expressed in HEK and CHO cells were not similar. This apparent discrepancy may be explained by the differences in the level of glycosylation of proBNP expressed in HEK and CHO cells.




FIGURE 2. Immunoreactivity of 6 external calibrators measured with 5 commercial BNP assays. The values of 4 measurements (0.585, 0.195, 0.065, 0.022 nM for each calibrator) were averaged (\pm SD).

TABLE 1. Agreement between BNP concentrations measured by different BNP assays with different calibrators.

Assays combination	Internal calibrator	Synthetic BNP	Recombinant BNP	Recombinant proBNP nonglyc	Recombinant proBNP nonglyc His-tagged	Recombinant proBNP glyc (CHO cells)	Recombinant proBNP glyc (HEK cells)
Alere Triage/Abbot i-STAT	$0.5238x + 3.5932$ $R^2 = 0.9897$	$0.7934x - 13.446$ $R^2 = 0.9886$	$0.8034x - 4.9589$ $R^2 = 0.9893$	$0.7934x - 13.446$ $R^2 = 0.9886$	$1.6224x + 33.221$ $R^2 = 0.9899$	$0.5462x - 6.3346$ $R^2 = 0.9893$	$1.007x - 33.151$ $R^2 = 0.9876$
Beckman Access2/Abbot i-STAT	$0.8063x - 53.328$ $R^2 = 0.9808$	$0.6991x - 37.912$ $R^2 = 0.97712$	$0.635x - 10.066$ $R^2 = 0.98383$	$0.6991x - 37.912$ $R^2 = 0.9771$	$0.9155x - 163.99$ $R^2 = 0.97263$	$0.6544x - 67.62$ $R^2 = 0.97605$	$1.0416x - 65.273$ $R^2 = 0.97769$
Beckman Access2/Alere Triage	$1.5334x - 56.349$ $R^2 = 0.98349$	$0.8792x - 25.309$ $R^2 = 0.98381$	$0.7873x - 5.2601$ $R^2 = 0.986867$	$0.6308x - 60.7$ $R^2 = 0.98162$	$0.5626x - 179.32$ $R^2 = 0.97656$	$1.1946x - 58.206$ $R^2 = 0.98288$	$1.032x - 29.689$ $R^2 = 0.98538$
Beckman Access2/Siemens Centaur XP	$1.366x - 25.375$ $R^2 = 0.99669$	$1.1755x - 22.718$ $R^2 = 0.99586$	$0.985x - 6.6407$ $R^2 = 0.99705$	$1.1682x - 21.734$ $R^2 = 0.99678$	$1.4711x - 46.968$ $R^2 = 0.99669$	$0.5641x - 21.809$ $R^2 = 0.99691$	$0.9567x - 63.68$ $R^2 = 0.99398$
ET Healthcare Pylon/Abbot i-STAT	$0.8237x + 33.157$ $R^2 = 0.95886$	$0.6396x + 43.999$ $R^2 = 0.9539$	$0.8656x + 21.854$ $R^2 = 0.95625$	$0.6991x - 37.912$ $R^2 = 0.9771$	$1.127x - 50.003$ $R^2 = 0.96059$	$0.5434x + 14.709$ $R^2 = 0.95895$	$0.927x + 26.061$ $R^2 = 0.95885$
ET Healthcare Pylon/Alere Triage	$1.5457x + 39.002$ $R^2 = 0.9361$	$0.789x + 61.157$ $R^2 = 0.92428$	$1.0565x + 33.313$ $R^2 = 0.92941$	$0.822x + 32.263$ $R^2 = 0.93539$	$0.6833x - 50.564$ $R^2 = 0.93883$	$0.9767x + 30.543$ $R^2 = 0.93418$	$0.9006x + 67.279$ $R^2 = 0.92915$
ET Healthcare Pylon/Beckman Access2	$0.999x + 101.23$ $R^2 = 0.93486$	$0.8868x + 87.051$ $R^2 = 0.91731$	$1.3381x + 41.219$ $R^2 = 0.93649$	$0.6991x - 37.912$ $R^2 = 0.9771$	$1.1955x + 184.97$ $R^2 = 0.93151$	$0.8096x + 82.674$ $R^2 = 0.93193$	$0.8665x + 96.399$ $R^2 = 0.92969$
ET Healthcare Pylon/Siemens Centaur XP	$1.3686x + 74.041$ $R^2 = 0.93725$	$1.0503x + 64.744$ $R^2 = 0.92743$	$1.32x + 31.873$ $R^2 = 0.93651$	$1.5095x + 88.746$ $R^2 = 0.93378$	$1.7643x + 125.15$ $R^2 = 0.93423$	$0.4579x + 63.812$ $R^2 = 0.93375$	$0.8369x + 36.367$ $R^2 = 0.94185$
Siemens Centaur XP/Abbot i-STAT	$0.5914x - 21.409$ $R^2 = 0.98794$	$0.5967x - 13.843$ $R^2 = 0.98749$	$0.6455x - 3.7951$ $R^2 = 0.98937$	$0.5967x - 13.843$ $R^2 = 0.9875$	$0.6244x - 82.068$ $R^2 = 0.98247$	$1.1627x - 83.832$ $R^2 = 0.98567$	$1.0923x - 3.6686$ $R^2 = 0.99005$
Siemens Centaur XP/Alere Triage	$1.1235x - 23.071$ $R^2 = 0.98838$	$0.7486x - 2.4625$ $R^2 = 0.98975$	$0.7994x + 1.3667$ $R^2 = 0.99007$	$0.5406x - 33.978$ $R^2 = 0.98737$	$0.3833x - 91.687$ $R^2 = 0.98426$	$2.1198x - 65.667$ $R^2 = 0.98797$	$1.0783x + 35.769$ $R^2 = 0.99053$

Almost 2-fold reduction in mean between-assay CV (%) was observed after recalibration with glycosylated proBNP expressed in HEK cells compared to internal calibrators (14.8% vs. 28.9%) (Fig. 3).




FIGURE 3. Equivalence of 5 commercial BNP assays calculated as the mean between-assay CV (%) for internal calibrators and 6 external calibrators (\pm SD).

Conclusions

1. Recombinant glycosylated proBNP (expressed in HEK cells) could serve as a common calibrator for BNP immunoassays to reduce between-assay variability and achieve better comparability of BNP concentrations of commercial BNP immunoassays.

2. This study may be considered as a starting point toward standardization of BNP measurement results. Further studies organized by standardization committees are needed to evaluate the potential of the suggested calibrator to become a reference material.