



## **Cardiac Markers**

Antibodies and antigens



# Introduction – cardiac biomarkers and diagnostics


**Cardiovascular diseases (CVDs) are a group of disorders of the heart and blood vessels. CVDs are the leading cause of death globally; it is estimated that approximately 30% of all deaths are caused by CVDs.**

First markers for cardiac diseases diagnostics were described already in the late 1950s and early 1960s when it was shown that measurements of lactate dehydrogenase (LDH), creatine kinase (CK) or aspartate aminotransferase (ASAT) enzymatic activities could be used in diagnosis of acute myocardial infarction (AMI). However, these assays had low specificity and sensitivity. First improvement was the development of immunological assays using polyclonal antibodies. In the early 1980s the first monoclonal antibodies were brought to the market marking a major leap in the development of immunoassays.

Since those days cardiac diseases diagnostics has gone through a significant evolution. Enzymatic assays and myoglobin assays have been replaced with cardiac troponins (I and T) and their high sensitivity versions are becoming the leading paradigm in AMI diagnostics. At the same time, markers for heart failure diagnosis have been adopted into routine use in most clinical laboratories. Within the recent years, more emphasis has been put on the development of markers which could be used for CVD prevention and risk assessment.


Hytest has been at the forefront of development of reagents for CVDs' diagnostics. During the past 20 years our scientists have authored or co-authored several articles published in peer reviewed scientific journals. This investment in scientific work has helped us to develop raw materials used by world's leading diagnostics companies. Hytest scientists have also worked actively in IFCC and AACC standardization committees. In 2004, Hytest's cardiac troponin complex material was chosen as a raw material for the international troponin standard.

We believe that reliable raw materials for diagnostics can be produced only if the design is based on solid understanding of the behavior of the biomarker, the disease state and most importantly, the needs of the industry using the raw materials.



Markers of myocardial infarction and heart failure. This schematic representation shows how they differ in timing and in specificity.

Introduction – cardiac biomarkers and diagnostics .....	2	C-reactive protein (CRP) .....	20
Cardiac troponin I (cTnI).....	4	Myoglobin .....	22
Cardiac troponin T (cTnT) .....	8	Myeloperoxidase (MPO).....	23
Human proBNP and its derivatives NT-proBNP and BNP.....	10	Fatty acid binding protein (FABP).....	24
ST2 .....	12	Other markers of cardiovascular disease...	25
Lipoprotein-associated phospholipase A2 (Lp-PLA2) .....	14	References.....	27
Pregnancy Associated Plasma Protein A (PAPP-A).....	16	Selected cardiac marker articles from Hytest scientists .....	28
Insulin-like growth factor binding protein-4 (IGFBP-4) and its fragments.....	18	Patents and trademarks.....	29
		More information.....	30



# Cardiac troponin I (cTnI)

## CLINICAL UTILITY

- Acute myocardial infarction (AMI)
- Unstable angina
- AMI prognosis
- Cardiac muscle injury and cell death

Cardiac troponin I is currently considered to be the gold standard biomarker test for myocardial infarction. Moreover, cTnI measurements by a new generation of high-sensitivity cTnI assays could be helpful for long-term risk stratification of different patient groups, including patients with heart failure or stable coronary artery disease.


At Hytest, we have intensively studied troponin I for over 30 years. Based on this research, we constantly aim to develop improved antibodies to be used in the immunoassays that are needed for accurate cardiac disease diagnostics. We have generated and tested several thousand monoclonal antibodies

specific to different regions of the cTnI molecule and have tested numerous different MAb combinations in order to find the best pairs for a precise and sensitive cTnI immunoassay.

## Factors influencing epitope recognition by antibodies


The most common reason for discrepancy in the cTnI assay measurements is the difference in the epitope specificity of the antibodies used in various assays. Due to several possible posttranslational modifications of the cTnI molecule found in patients' blood and the presence of autoantibodies in some clinical samples, it is critical to carefully validate the performance of antibodies in order to achieve reliable, quantitative detection of cTnI in blood samples.

When designing a sensitive and precise immunoassay, it is important to consider the effect of all the factors influencing biomarker detection. The assay should not be affected by partial proteolytic degradation of cTnI molecule, oxidation, reduction, phosphorylation, complex formation with TnC, or the presence




**Figure 1.** Epitope mapping of Hytest anti-cTnI monoclonal antibodies. We offer several specially selected antibodies specific to various epitopes along the cTnI molecule.


# At Hytest, we have intensively studied troponin I for over 30 years.



**Figure 2.**  
Factors influencing cTnI immunodetection.



**Figure 3.**  
**Effect of proteolytic degradation.** Best two-site combinations of cTnI antibodies specific to the stable part of cTnI molecule tested with troponin complex before (light red columns) and after (blue columns) incubation for 170 hours with a mixture of endogenous proteases from human cardiac tissue. Control assay M18-MF4 is sensitive to cTnI proteolytic degradation.




**Figure 4.**  
**Effect of heparin.** cTnI concentration was measured in the absence (light red columns) or presence (blue columns) of 5 IU/ml heparin. Antibody 228 in the assay 228-19C7 is sensitive to the presence of heparin in the sample.

of heparin in the samples (Katrukha, 2003). This helps to minimize bias in the assay. Factors that influence cTnI measurements are schematically presented in Figure 2 and examples of the effects of these interfering factors on analyte detection are shown in Figures 3 and 4.

Antibodies specific to different parts of the molecule are sensitive to interfering factors in different degrees. For instance, it is well known that purified cTnI is highly susceptible to proteolytic degradation. However, in a troponin complex the central part of the cTnI closely interacts with TnC which protects cTnI from endogenous proteases. Consequently, the epitopes located in the central part of the cTnI are significantly more stable than the epitopes located at the terminal parts of the molecule. On the other hand, not every antibody specific to the central part of the molecule can recognize cTnI in a patient's blood because TnC covers some of the epitopes located in that region.

## Antibodies for high-sensitivity cTnI immunoassays

In an immunoassay, the limit of detection is dependent on many features — platform, label, incubation time, buffers used, and many others. However, the most critical is the affinity of antibodies that are used for the assay design. Today Hytest antibodies are successfully utilized for the development of a new generation of high-sensitivity cTnI assays. An example of a highly sensitive assay is shown in Figure 5.



**Figure 5.**  
**Highly sensitive quantitation of cTnI (down to 0.64 pg/ml) using MAb 801 and 19C7.** With LamdaGen's plasmonic ELISA platform utilizing MABs 801 and 19C7 from Hytest. The limit of detection for cTnI was 0.64 pg/ml. The dose response curve was obtained by spiking native cTn complex into 87% fetal bovine serum. Each data point represents the average of five independent measurements. Reprinted with permission from LamdaGen Corp.  
<sup>1</sup>OES™ quantification of cTnI in whole serum. Application Note, 2013. LamdaGen Corporation. www.LamdaGen.com

## Chimeric anti-cTnI antibodies

Heterophile antibodies arise when people are exposed to different animals or products derived from animals. As far as immunodiagnostics is concerned, the problem is most commonly associated with human anti-mouse antibodies (HAMA) due to the fact that most diagnostics assays use mouse derived antibodies. HAMA might cause both false negative and false positive results that could lead to delays in making the correct diagnosis. Troponin assays are particularly susceptible to HAMA due to low cut-off value requirements and because the levels of cTnI even in the plasma of AMI patients are very low.

A powerful tool to solve the issue with HAMA in diagnostics tests is the use of chimeric or fully humanized antibodies. We have now converted two of our cTnI antibodies, MAbs 19C7 and 16A11, to chimeric proteins by changing the antibody constant regions from mouse to human derived sequences. The chimeric cTnI antibodies RecChim19C7 and RecChim16A11 consist of the original mouse derived variable regions that are responsible for antigen specificity and human derived constant regions of IgG1 isotype (see Figure 6).

## Chimeric antibodies prevent the HAMA effect

The performance of different combinations of chimeric and native antibodies was tested using HAMA containing serum samples that were obtained from acute myocardial patients in order to verify that the chimeric antibodies are not sensitive to the HAMA effect (see Figure 7).

## Hytest troponin complex selected as reference material

cTnI, which is extremely unstable in its free form, demonstrates significantly better stability in complex with TnC or in ternary cTnI-cTnT-TnC complex (not shown). These two forms of the protein are preferable as material for standard and calibrator preparation. In the troponin complex supplied by Hytest, cTnI is presented in the same form as it can be detected in the blood of AMI patients. Purification of the troponin complex is performed in mild conditions without treatment with urea containing buffers (as is usually done when preparing individual troponin components). The concentration is precisely determined for each of the three components of the complex.




Figure 6. Schematic illustration of Hytest chimeric cTnI antibodies.




Figure 7. Chimeric antibodies mitigate the HAMA effect. The performance of chimeric and native 19C7 and 16A11 in the presence of HAMA was tested with three serum samples with varying HAMA concentrations: 807 ng/ml in Sample 1, 1388 ng/ml in Sample 2 and 6220 ng/ml in Sample 3. As a control, buffer without serum was used. Antibody pairs compared are indicated in the picture.

Hytest troponin complex selected as reference material by various standardization bodies.

## MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4T21cc* 4T21*	Monoclonal mouse anti-cardiac troponin I (cTnI)	Enzyme immunoassays Western blotting Immunoprecipitation Immunohistochemistry Immunoaffinity purification
RC4T21*	Recombinant anti-cardiac troponin I (cTnI)	Enzyme immunoassays
4T45	Monoclonal mouse anti-cardiac troponin I (cTnI), phosphorylated form	Enzyme immunoassays Western blotting
4T46	Monoclonal mouse anti-cardiac troponin I (cTnI), dephosphorylated form	Enzyme immunoassays Western blotting
4TC2*	Monoclonal mouse anti-cardiac troponin complex	Enzyme immunoassays
RC4TC2	Recombinant anti-cardiac troponin complex	Enzyme immunoassays
4T20*	Monoclonal mouse anti-skeletal muscle troponin I (skTnI)	Enzyme immunoassays Western blotting

\* Several MAbs available under one catalogue number. Please see [www.hytest.fi](http://www.hytest.fi).

## POLYCLONAL ANTIBODY

Cat.#	Product	Host	Tested applications
4T21/2*	Polyclonal anti-cardiac troponin I (cTnI)	Goat	Immunoassays

\* Available in certain regions.

## ANTIGENS

Cat.#	Product	Source	Purity
8T53	Troponin I cardiac, human	Human cardiac muscle	>98%
8RT17	Troponin I cardiac, human, recombinant	Recombinant	>95%
8T25	Human skeletal TnI	Human skeletal muscle	>95%
8ICR3	Troponin IC complex, cardiac, human	Recombinant	>95%
8IFC20	Troponin I (fragment 28-110) – troponin C complex, cardiac, human, recombinant chimeric	Recombinant	>95%
8ITCR	Troponin ITC complex, cardiac, human	Recombinant	>95%

## SERUM AND OTHER PRODUCTS

Cat.#	Product	Source/Remarks
8TFS	cTnI free serum	Pooled normal human serum

# Cardiac troponin T (cTnT)

## CLINICAL UTILITY

- AMI
- Unstable angina
- AMI prognosis
- Cardiac muscle injury and cell death

Cardiac troponin T (cTnT), along with cardiac troponin I (cTnI) is accepted as a “golden marker” for myocardial infarction (MI) diagnosis. Both biomarkers are released into circulation with same kinetics and either of them can be used for the diagnosis of MI when tested using contemporary or point of care instruments. Moreover, cTnT and cTnI measurements by high-sensitivity assays could be helpful for long-term risk stratification of different patient groups with cardiac disease and/or for early rule-out or rule-in patients in emergency departments.

## High-sensitivity cTnT assay prototypes

Assay prototypes utilizing our newly developed anti-cTnT monoclonal antibodies demonstrate a good linearity and superior sensitivity. The limit of detection (LoD) of both assays was better than 0.3 ng/l. A typical calibration curve for a purified cTnT in low concentration (0.14-100 ng/l) is presented in Figure 8.

## MI blood testing and correlation with a commercially available hs-cTnT assay

We studied the correlation of the two prototype immunoassays to a commercially available hs-cTnT assay by analyzing 38 serum samples from AMI patients. Figure 9 shows that there is a good correlation between cTnT values obtained with our assay prototypes and the commercially available hs-cTnT assay. The new cTnT MAbs allow for the development of highly sensitive immunoassays for the detection of cTnT in the blood of AMI patients with high specificity.




Figure 8.

Calibration curve for a MAb combination 329-406. Purified native human cTnT was used as the antigen.




Figure 9.

Hytest immunoassays show good correlation to a commercially available hs-cTnT assay. Concentration of cTnT in 38 serum samples obtained from AMI patients was determined by using two immunoassays utilizing Hytest antibodies (capture-detection pairs 329-406 and 406-300) and a commercially available hs-cTnT assay.

## MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4T19cc* 4T19*§	Monoclonal mouse anti-cardiac troponin T (cTnT)	Western blotting Affinity purification Immunohistochemistry Immunoprecipitation
RC4T19	Recombinant anti-cardiac troponin T (cTnT)	Enzyme immunoassays

\* Several MAbs available under one catalogue number. Please see [www.hytest.fi](http://www.hytest.fi).

§ Recommended for research purposes.

## ANTIGENS

Cat.#	Product	Source	Purity
8RTT5	Human cardiac TnT, recombinant	Recombinant	>95%
8T24	Human skeletal TnT	Human skeletal muscle	>95%
8RTF4	Recombinant human fast skeletal troponin T (fast skTnT)	Recombinant	>95%
8RST2	Recombinant human slow skeletal troponin T (slow skTnT)	Recombinant	>95%

# Human proBNP and its derivatives NT-proBNP and BNP

## CLINICAL UTILITY

- Identification or exclusion of heart failure (HF)
- Assessment of the severity of HF
- Prognosis of the disease development
- Monitoring of drug therapy in the presence of HF

Pro-B-type natriuretic peptide (proBNP) derivatives BNP and NT-proBNP are established biomarkers in heart failure (HF) diagnostics (Figure 10). The concentrations of BNP and NT-proBNP in blood increase rapidly as a consequence of cardiac wall stretch and correlate with the severity of the disease. Analysis of BNP and NT-proBNP levels are used e.g. for exclusion of heart failure, risk stratification and as a prognostic marker of heart failure.

Information obtained from BNP and NT-proBNP studies conducted within the last few years has greatly improved our understanding of the properties and processing of proBNP. This, in turn, makes it easier to design immunoassays which


reliably and quantitatively detect the biomarkers from clinical samples.

## NT-proBNP detection is affected by glycosylation


Our studies have revealed that the majority of antibodies specific to the central part of the NT-proBNP molecule scarcely detect the antigen in human blood samples. We found that this is due to glycosylation of the central part of endogenous NT-proBNP. Glycosylation changes the epitope availability for antibody recognition, subsequently decreasing reliable quantitation of NT-proBNP with such antibodies (Figure 11 and Seferian et al., 2008). For precise NT-proBNP measurements in human blood, we recommend using a pair of antibodies specific to the N- and C-terminal parts of the NT-proBNP molecule.

## BNP is an unstable molecule

BNP belongs to a family of peptide hormones that all have a 17 amino acid ring structure with a disulfide bond between two cysteine residues. BNP is an unstable molecule and it contains several protease cleavage sites. Especially the N-terminus is highly susceptible to degradation. The ring structure is relatively stable. In order to reliably detect BNP, it is recommended to choose antibodies that are less sensitive to proteolytic degradation of BNP.



**Figure 10.** Schematic representation of proBNP processing. ProBNP is processed in a convertase-dependent reaction into NT-proBNP and BNP. BNP has biological activity whereas the role of NT-proBNP is unknown.



**Figure 11.**

**Immunoreactivity of endogenous NT-proBNP before and after deglycosylation.** Concentration of endogenous NT-proBNP before (light red columns) and after (blue columns) deglycosylation was measured in 8 HF patient samples by sandwich immunoassays. (A) MAbs 15C4<sub>63-71</sub> and 13G12<sub>13-24</sub> are specific to the N- or C-terminal parts of the molecule which are not glycosylated. (B) The capture MAb (11D1<sub>31-39</sub>), recognizes an epitope located at the central region. This region of endogenous NT-proBNP becomes available for antibody recognition only after deglycosylation.

## Hytest NT-proBNP selected as reference material by The National Institute of Metrology (NIM) of China

### MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4BNP2cc* 4BNP2	Monoclonal mouse anti-human brain natriuretic peptide (BNP)	Enzyme immunoassays Western blotting
4BFab5cc	Monoclonal mouse anti-immune complex (24C5-BNP/proBNP)	Enzyme immunoassays
4NT1cc* 4NT1*	Monoclonal mouse anti-human N-terminal proBNP (NT-proBNP)	Enzyme immunoassays Western blotting

\* Several MAbs available under one catalogue number. Please see [www.hytest.fi](http://www.hytest.fi).

### ANTIGEN

Cat.#	Product	Source	Purity
8NT2	NT-proBNP, recombinant	Recombinant	>95%
8PRO9	ProBNP, recombinant	Recombinant	>95%
8GBP3	ProBNP, glycosylated, recombinant	Recombinant	>95%

### DEPLETED PLASMA

Cat.#	Product	Source/Remarks
8BFP	BNP and NT-proBNP free plasma	Pooled normal human plasma

# ST2 – a marker of cardiac stress

## CLINICAL UTILITY

- **Prognostic marker of heart failure**

ST2 exists in two isoforms: a transmembrane or cellular (ST2L) and soluble or circulating (sST2). The IL-33 exerts its effects by binding to the transmembrane receptor ST2L isoform, and the IL-33/ST2 system is upregulated in cardiomyocytes and fibroblasts in response to cardiac injury. The expression of sST2 is largely inducible and it is almost ubiquitous in living cells, such as resting fibroblasts. It has been suggested that sST2 is produced by both cardiac fibroblasts and cardiomyocytes in response to cardiac injury or cardiac stress, macrovascular (aortic and coronary artery), and heart microvascular endothelial cells in response to diastolic load.

## Clinical value of ST2

A sST2 is a biomarker that is used for additive risk stratification and prognosis of patients with heart failure (HF). In contrast to BNP and NT-proBNP, ST2 is not affected by confounding factors such as age, body mass index, and impaired renal function. Unlike many other cardiac biomarkers, the levels of ST2 alter quickly in response to changes in the patient's condition. Elevated levels of sST2 in acute as well as chronic HF patients (>35 ng/ml) are strongly associated with the measures of HF severity, and they predict both recurrent hospitalization and mortality.

Measured levels of ST2 in chronic HF patients can be used for therapy evaluation and accordingly, decreased levels of sST2 that are responsive to medical treatment are associated with better outcomes for patients. According to most publications, the ST2 is an independent predictor of all-cause cardiac mortality, and it provides complementary prognostic information not only for NT-proBNP (or BNP), but also for high-sensitivity cardiac troponin T (hs-cTnT) assays.


## Reagents for the ST2 immunoassay development

Hyttest ST2 MABs have been studied in different immunoassay formats, such as sandwich ELISA and immunofluorescent (FIA) assays. Both assay platforms can be used for the development of high-sensitive and precise ST2-specific immunoassays, with a level of detection (LoD) up to 30 pg/ml.


Hyttest recombinant ST2 expressed in mammalian cell line is suitable for use as a standard or calibrator in immunoassays. The protein has a theoretical molecular weight (Mw) of approximately 36 kDa and it has several sites of N-glycosylation.

All recommended capture-detection pairs that we have tested in the sandwich ELISA platform have sensitivity that is comparable with the sensitivity of the reference Presage® ST2 assay or better, while using serial dilutions of pooled plasma from HF patients. A calibration curve for the prototype assay S215-S103 is shown in Figure 12.

Our capture-detection pairs have been tested in the sandwich immunoassay and the sensitivity is comparable with the sensitivity of the reference Presage® ST2 assay or better. Figure 13 shows the correlation studies between immunoassays that utilize MABs S215-S103 and a commercially available diagnostic Presage® ST2 assay. The concentrations of ST2 obtained with these assays are directly comparable with the Pearson correlation coefficient of 0.98.



**Figure 12.**  
**Representative calibration curve for the ST2 prototype sandwich ELISA assay (S215-S103), using recombinant ST2 as the antigen.**



**Figure 13.**  
**ST2 concentrations that were obtained with the S215-S103 prototype immunoassay and the Presage® ST2 assay.** The correlation coefficient (Pearson) between the assay S215 (capture) - S103 (detection) and the Presage ST2 assay is 0.98. Both assays are sandwich ELISA immunoassays and the protocol used was similar to that which was used in the Presage® ST2 assay.

## MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4ST2*	Anti-ST2	Enzyme immunoassays FIA

\* Several MAbs available under one catalogue number. Please see [www.hytest.fi](http://www.hytest.fi).

## ANTIGEN

Cat.#	Product	Source	Purity
8STR4	ST2/ IL1RL1 protein, human, recombinant	Recombinant	>95%

# Lipoprotein-associated phospholipase A2 (Lp-PLA2)

## CLINICAL UTILITY

- **Prognostic marker of adverse cardiac-related events**

Lipoprotein-associated phospholipase A2 (Lp-PLA2, which is also known as the platelet-activating factor acetyl hydrolase) is a  $\text{Ca}^{2+}$ -independent phospholipase that circulates in the bloodstream in the form of a complex with lipoprotein particles (Stafforini, 2009; 2015). Lp-PLA2 levels have been shown to predict adverse cardiac-related events in both patients with stable coronary artery disease (Brilakis et al., 2005) and in a healthy adult population (Ballantyne et al., 2004). The increase in Lp-PLA2 levels can predict the development of incident peripheral arterial disease in humans (Garg et al., 2016).

Recent guidelines from four major international societies, which include the European Society of Cardiology, the American College of Cardiology, the American Heart Association and the American Society of Endocrinology, have included Lp-PLA2 among the biomarkers whose measurement could be useful for risk stratification of asymptomatic adult patients.


## Recombinant human Lp-PLA2

Hyttest provides recombinant human Lp-PLA2 (recLp-PLA2) that is expressed in a mammalian cell line. The protein contains 6×His tag on its C-terminus linked with a GG spacer.

Upon serial dilutions, recombinant Lp-PLA2 and endogenous Lp-PLA2 in normal human serum showed the same pattern of signal decrease in sandwich fluoroimmunoassays employing the MAb combination PL42cc-PL46cc (see Figure 14). This demonstrates that the immunochemical properties of recombinant Lp-PLA2 are similar to those of native Lp-PLA2.

## Measuring patient samples

We provide several different monoclonal antibodies specific to human Lp-PLA2 that allow the development of quantitative immunoassays for detecting endogenous Lp-PLA2 in serum samples. When compared to a commercially available ELISA kit, the Hyttest antibodies detected Lp-PLA2 in a manner that was very similar or slightly different to the commercial kit (see Figure 15).




**Figure 14.**

**Dilutional linearity study.** Dilutional linearity study of recombinant Lp-PLA2 and native Lp-PLA2 (normal human serum from an apparently healthy volunteer) studied using the MAb combination PL42cc-PL46cc. The initial concentration of the recombinant human Lp-PLA2 was 111 ng/ml.

**Figure 15.**

**Detection of native Lp-PLA2 in serum samples.** Lp-PLA2 was detected in serum of acute myocardial infarction patients and healthy volunteers using the Hytest fluoroimmunoassay with the MAb pair PL26cc-PL4cc (A) and PL42cc-PL46cc (B) or by using a commercially available ELISA kit (C). Serum samples were diluted 1:30 with an assay buffer (A and B) or according to the manufacturer's instructions (C). Samples were incubated for 2.5 hours at 37°C (A and B) or for 3 hours at room temperature (C).



## MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4LA7cc*	Monoclonal mouse anti-human lipoprotein-associated phospholipase A2 (Lp-PLA2)	Enzyme immunoassays Western blotting

\* Several MABs available under one catalogue number. Please see [www.hytest.fi](http://www.hytest.fi).

## ANTIGEN

Cat.#	Product	Source	Purity
8PL7	Recombinant human lipoprotein-associated phospholipase A2 (Lp-PLA2)	Recombinant	>75%

# Pregnancy Associated Plasma Protein A (PAPP-A)

## CLINICAL UTILITY

- Acute myocardial infarction
- Acute coronary syndrome
- Unstable angina
- Down syndrome

The pregnancy-associated plasma protein-A (PAPP-A) has been used as a biochemical marker for Down syndrome in the first trimester of pregnancy for a long time. In addition to this, several studies show that PAPP-A is a promising marker for cardiac diseases.

PAPP-A can be found in blood in two different forms: as a heterotetrameric complex (htPAPP-A) and as a homodimeric complex (dPAPP-A) (Figure 16). Of these, the dimeric form is the one associated with cardiac diseases. dPAPP-A is abundantly expressed in unstable coronary atherosclerotic plaques (Bayes-Genis et al., 2001).

It has been demonstrated that the level of dPAPP-A in the blood is significantly elevated in patients with unstable angina or acute myocardial infarction, in comparison to patients with stable angina and control subjects (Heeschen et al., 2005, Hájek et al., 2008). In addition, dPAPP-A has also been shown to be a strong independent marker of risk stratification for patients with acute coronary syndrome (ACS) (Qin et al., 2002). In their literature search report, Richard Body and Craig Ferguson (2006) concluded that PAPP-A is a promising biomarker for unstable coronary disease and that it also could have great potential as a prognostic marker as part of a multimarker strategy.

## dPAPP-A specific sandwich immunoassay

Our dPAPP-A specific monoclonal antibodies enable the development of immunoassays that are suitable for selective quantitative measurements of dPAPP-A in human blood, even in the presence of htPAPP-A (Figure 17).




Figure 16.

**Structure of PAPP-A.** PAPP-A is a metalloprotease that belongs to the metzincin superfamily of zinc peptidases. Homodimeric dPAPP-A consists of two 200 kDa PAPP-A subunits covalently linked with a disulfide bond. Heterotetrameric htPAPP-A includes two PAPP-A subunits and two 50-90 kDa subunits of the preform of the eosinophil major basic protein (proMBP), all covalently linked with disulfide bonds. It has been shown that proMBP has inhibitory properties against the protease activity of PAPP-A.




Figure 17.


**dPAPP-A measured in clinical samples.** Concentration of dPAPP-A in plasma samples of 43 ACS patients (ACS) and 34 non-ACS patients control group (Normal) measured by PAPP52 – PAPP30 sandwich immunoassay (mean $\pm$ SD).

Capture MAb: PAPP52 (Cat.# 4P41) Detection MAb: PAPP30 (Cat.# 4PD4; labeled with Eu<sup>3+</sup> chelate)

Incubation volume: 100  $\mu$ l. Incubation time: 30 min at room temperature.

## Dimeric recombinant human PAPP-A

Our dimeric recombinant human PAPP-A contains His10-tag and is expressed in mammalian cells and purified by metal affinity chromatography (Figure 18). The product has E483A mutation for stabilization of the protein due to suppression of proteolytic activity and autocleavage.



**Figure 18.**

### SDS-PAGE of dimeric recombinant human PAPP-A.

1. in non-reducing conditions, 5 µg;
2. molecular weight markers;
3. in reducing conditions, 5 µg.

The level of dPAPP-A in the blood is significantly elevated in patients with unstable angina or acute myocardial infarction.

## MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4P41*	Monoclonal mouse anti-human pregnancy-associated plasma protein A (PAPP-A)	Enzyme immunoassays Western blotting
4PD4	Monoclonal mouse anti-human dimeric form of pregnancy-associated plasma protein A (dPAPP-A)	Enzyme immunoassays

\* Several MAbs available under one catalogue number. Please see [www.hytest.fi](http://www.hytest.fi).

## ANTIGEN

Cat.#	Product	Source	Purity
8PA1	PAPP-A, human, recombinant	Recombinant	>90%

# Insulin-like growth factor binding protein-4 (IGFBP-4) and its fragments

## CLINICAL UTILITY

- Prediction of major adverse cardiac events
- Acute myocardial infarction
- Acute coronary syndrome
- Unstable angina

IGFBP-4 has been shown to be a substrate for dPAPP-A (Figure 19). dPAPP-A is a promising marker for predicting plaque rupture which, in turn, may lead to acute thrombosis. However, measuring dPAPP-A concentrations reliably is challenging due to many reasons (Terkelsen et al., 2009, Tertti et al., 2009). Our studies indicate that quantifying N- and C-terminal IGFBP-4

fragments (NT-IGFBP-4 and CT-IGFBP-4, respectively) that are the result of dPAPP-A cleavage instead of dPAPP-A could be used as an indirect but more reliable method for obtaining information about dPAPP-A concentration and, consequently, for predicting the rupture of vulnerable plaques (Postnikov et al., 2012).





Figure 19.

Schematic representation of dPAPP-A dependent cleavage of IGFBP-4 in unstable atherosclerotic plaque. dPAPP-A cleaves IGFBP-4 (preferably complexed with IGF) into NT-IGFBP4 and CT-IGFBP-4 fragments. As a result, IGF is released and activated. IGFBP-4 fragments are released into circulation and can be detected from blood.

## Immunoassays for quantifying NT-IGFBP4 and CT-IGFBP4

In order to quantify IGFBP-4 fragments, we have generated MAbs specific to epitopes available for MAb binding only after the proteolytic cleavage of IGFBP-4 by dPAPP-A (Figure 20). Cross-reactivity of these neoepitope-specific MAbs with full-length IGFBP-4 is negligible (1.4% or less), thus allowing for specific quantitation of cleaved fragments regardless of the presence of non-cleaved IGFBP-4.



**Figure 20.** Sandwich immunoassays for detecting NT- and CT-IGFBP-4. (A) Schematic representation of the assays and capture/detection MAbs chosen for each assay. (B) Representative immunoassays with purified recombinant NT- and CT-IGFBP-4 fragments.

We have generated MAbs specific to epitopes available for MAb binding only after the proteolytic cleavage of IGFBP-4 by dPAPP-A.

### MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4IGF4*	Monoclonal mouse anti-Insulin-like growth factor binding protein 4 (IGFBP-4)	Enzyme immunoassays

\* Several MAbs available under one catalogue number. Please see [www.hytest.fi](http://www.hytest.fi).

### ANTIGEN

Cat.#	Product	Source	Purity
8NGP4	Insulin-like growth factor binding protein 4, N-terminal fragment (NT-IGFBP-4), human	Recombinant	≥90%
8ILG4	Insulin-like growth factor binding protein 4, C-terminal fragment (CT-IGFBP-4), human	Recombinant	≥90%

# C-reactive protein (CRP)

## CLINICAL UTILITY

- **Prediction of future cardiovascular risk**
- **Inflammation**

C-reactive protein (CRP) is one of the so-called acute phase proteins. Its concentration in blood increases rapidly as a response to inflammation.

In recent years, more information has been obtained about the possible role of inflammation in contributing to the development of serious health issues such as diabetes or the development of cardiovascular diseases. These studies show that elevated basal levels of CRP indicate increased risk for cardiac diseases, thus making CRP a promising biomarker for predicting future development of a heart disease. Many epidemiologic studies have indicated that CRP is a strong independent predictor of future cardiovascular events, including myocardial infarction, ischemic stroke, peripheral vascular disease, and sudden cardiac death without a known cardiovascular disease (as reviewed by Clearfield, 2005).

## High sensitivity CRP (hsCRP) immunoassays

In 2003, the Centers for Disease Control and Prevention (CDC) and the American Heart Association (AHA) issued a statement that identified CRP as the inflammatory marker best suited for use in current clinical practice to assess cardiovascular risk (Ridker, 2003). While the CRP level in blood can rapidly increase to tens or even hundreds of milligrams per liter during an acute inflammation, it is the basal level of blood CRP that has more clinical significance when predicting future cardiac diseases (Scirica et al., 2007; Koenig et al., 2008). This is why present day hsCRP assays are aimed at nanogram per milliliter (ng/ml) CRP level distinction (Figure 21).





Figure 21.

**Immunodetection of CRP standard in a sandwich immunoassay by MAb pair C2-C6.** MAb C2 is biotinylated, MAb C6 is labeled with stable  $\text{Eu}^{3+}$  chelate. The mixture of antibodies and antigen samples (100  $\mu\text{l}$ ) was incubated for 10 min at room temperature in streptavidin coated plates.

## Antibodies for hsCRP immunoassays

Anti-CRP antibodies developed by Hytest have been utilized in several immunoassays achieving excellent sensitivity and a broad linear detection range (Meyer et al., 2007; Shiesh et al., 2006; Sin et al., 2006). These antibody combinations could be used for the development of hsCRP assays for different diagnostic platforms. In addition, for the convenience of our customers, we have monoclonal antibodies with different affinity (Figure 22 and Table 1), thus enabling them to be used in different types of immunoassays.




**Figure 22.**  
**Biacore X sensograms of different anti-CRP MAbs.** 100 nM native CRP was exposed to the chip-immobilized MAbs in HBS-EP buffer (0.01 M HEPES, 0.15 M NaCl, 3 mM EDTA, 0.005% polysorbate 20, pH 7.4).

**Table 1.**  
**Affinity constants for selected Hytest anti-CRP MAbs.**

MAb	$K_{on}$ (1/Ms)	$K_{off}$ (1/s)	$K_d$ (M)
C2cc	$2.3 \times 10^5$	$4.4 \times 10^{-4}$	$1.93 \times 10^{-9}$
C5	$1.3 \times 10^5$	$2.2 \times 10^{-3}$	$1.7 \times 10^{-8}$
CRP30cc	$9.3 \times 10^3$	$4.0 \times 10^{-3}$	$4.3 \times 10^{-7}$
CRP135cc	$1.5 \times 10^5$	$6.6 \times 10^{-4}$	$4.4 \times 10^{-9}$

## Recombinant human CRP

Hytest recombinant human CRP is expressed in mammalian cells and purified in native conditions that excludes renaturation steps. Recombinant human CRP is purified by affinity chromatography with phosphatidylcholine matrix that confirms functional activity of the recombinant protein. Our recombinant human CRP (Cat.# 8CR8) does not contain any tags. Recombinant human CRP is immunochemically active in different sandwich immunoassay pairs using Hytest antibodies (Figure 23).



**Figure 23.**  
**Calibration curve for the human recombinant CRP.** Calibration curve obtained with CRP30cc-CRP135cc (capture-detection) sandwich immunoassay.

## MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4C28cc* 4C28*	Monoclonal mouse anti-human C-reactive protein (high sensitivity CRP)	Enzyme immunoassays (also high sensitivity) Western blotting Turbidimetric assays Immunohistochemistry Immunoaffinity purification

\* Several MAbs available under one catalogue number. Please see [www.hytest.fi](http://www.hytest.fi).

## ANTIGEN

Cat.#	Product	Source	Purity
8CR8	C-reactive protein (CRP), human, recombinant	Recombinant	>95%

## DEPLETED SERUM


Cat.#	Product	Source/Remarks
8CFS	C-reactive protein free serum	Pooled normal human serum

# Myoglobin

## CLINICAL UTILITY

- Myocardial damage
- Acute myocardial infarction

Myoglobin has been used as a marker of myocardial damage for almost six decades (Figure 24). It is commonly used in clinical practice as an early marker of AMI (Penttilä *et al.*, 2002). However, due to the high concentration of myoglobin in skeletal muscle tissues, even minor skeletal muscle injury increases the myoglobin levels in blood (van Nieuwenhoven *et al.*, 1995). Therefore, myoglobin alone is not considered to be a reliable and sufficient marker in AMI diagnosis. Instead, it should be used together with cTnI or cTnT analysis as part of a multi-marker strategy.



**Figure 24.**  
SDS-PAGE of human myoglobin. M is molecular weight standard (Pharmacia).

## MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4M23*	Monoclonal mouse anti-human cardiac myoglobin	Enzyme immunoassays

\* Several MAbs available under one catalogue number. Please see [www.hytest.fi](http://www.hytest.fi).

## ANTIGEN

Cat.#	Product	Source	Purity
8M50	Myoglobin	Human cardiac muscle	>95%

## DEPLETED SERUM

Cat.#	Product	Source/Remarks
8MFS	Myoglobin free serum	Pooled normal human serum

# Myeloperoxidase (MPO)

## CLINICAL UTILITY


- Acute coronary syndrome
- Coronary artery disease
- Cardiovascular disease risk stratification
- Prediction of long term incident major adverse cardiac events

Myeloperoxidase (MPO) is a peroxidase enzyme that is abundantly secreted by activated leukocytes (neutrophils) during an inflammation reaction. Within the last decades, multiple studies have indicated that MPO is a promising cardiac marker. Brennan *et al.* (2003) showed that unlike troponins, CK-MB and CRP, MPO facilitates the identification of patients that are at risk for cardiac events in the absence of myocardial necrosis. It was also demonstrated that an increased level of MPO in a patient's blood serves as a risk marker for atherosclerosis (Nambi, 2005) and coronary artery disease (Zhang *et al.*, 2001).

Our anti-MPO antibodies have been tested with clinical samples.

## Sandwich immunoassay for quantitative MPO detection

All our MABs have been screened to provide sensitive and specific detection of endogenous MPO with good kinetics (Figure 25). In addition, we have tested our most sensitive two site MAB combinations with blood samples containing high titers of MPO autoantibodies and we have selected MAB combinations that are less sensitive to autoantibodies.



**Figure 25.** Calibration curve for a MPO sandwich immunoassay. MABs 16E3 and 18B7 were used as capture and detection antibodies, respectively. Native purified human MPO was used as the antigen.

## MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4M43*	Monoclonal mouse anti-human myeloperoxidase (MPO)	Enzyme immunoassays Western blotting

\* Several MABs available under one catalogue number. Please see [www.hytest.fi](http://www.hytest.fi).

## DEPLETED SERUM

Cat.#	Product	Source/Remarks
8MPFS	Myeloperoxidase free serum	Pooled normal human serum


# Fatty acid binding protein (FABP)

Fatty acid-binding proteins (FABPs) are a group of small cytoplasmic proteins that are abundant in tissues with active fatty acid metabolism. This includes the heart (Storch and Thumser, 2000). The heart-type fatty acid binding protein (H-FABP) is an early marker of myocardial injury and is widely

applied in emergency triage of patients with acute coronary syndromes (Alhadi and Fox, 2004). H-FABP is considerably more cardio-specific than myoglobin, another early AMI marker. However, it is less cardio-specific than troponins due to the fact that some H-FABP is also expressed in skeletal muscle tissues.

## Various MAb pairs can be used for H-FABP measurement from blood samples of AMI patients

Our anti-FABP antibodies allow for the development of quantitative, highly sensitive immunoassays for the detection of H-FABP (Figure 26). All MAbs have been tested with clinical samples.



**Figure 26.** Detection of H-FABP in serum samples of AMI patients. Comparison of H-FABP measurements in sera of six AMI patients using different antibody combinations in sandwich immunoassays.

## CLINICAL UTILITY

- Acute coronary syndrome
- Myocardial injury

## MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4F29*	Monoclonal mouse anti-human fatty acid binding protein (FABP)	Enzyme immunoassays Western blotting

\* Several MAbs available under one catalogue number. Please see [www.hytest.fi](http://www.hytest.fi).

## ANTIGEN

Cat.#	Product	Source	Purity
8F65	Fatty acid binding protein	Human cardiac muscle	>95%

## DEPLETED SERUM


Cat.#	Product	Source/Remarks
8FFS	Fatty acid binding protein free serum	Pooled normal human serum

# Other markers of cardiovascular disease

## Troponin C (TnC)

Troponin C (TnC) is the Ca<sup>2+</sup>-binding subunit of the troponin complex. In human muscle cells, it exists in two different isoforms, fast and slow. In myocardium, TnC is presented by the slow skeletal isoform.

TnC forms high affinity complexes with cTnI. It was demonstrated that cTnI is presented mainly as a complex with TnC in the blood stream of AMI patients (Katrukha *et al.*, 1997). In the binary cTnI-TnC complex, TnC protects cTnI from protease cleavage. Therefore, TnC can be used as a natural stabilizer of cTnI in water solutions (Katrukha *et al.*, 1998).



**Figure 27.** Native TnC. 1 µg of human native TnC (Cat.# 8T57) was run in 10-20% SDS-PAGE under reducing conditions.

### MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4T27cc	Monoclonal mouse anti-troponin C (TnC), <i>in vitro</i>	Enzyme immunoassays Western blotting
4TC2*	Monoclonal mouse anti-human cardiac troponin complex	Enzyme immunoassays
RC4TC2	Recombinant anti-cardiac troponin complex	Enzyme immunoassays

\* Several MAbs available under one catalogue number. Please see [www.hytest.fi](http://www.hytest.fi).

### ANTIGENS

Cat.#	Product	Source	Purity
8T57	Troponin C, human	Human cardiac muscle	>98%
8RSC4	Recombinant human slow skeletal/cardiac troponin C (TnC)	Recombinant	>95%
8RKC3	Recombinant human troponin C skeletal muscle, isoform 2	Recombinant	>90%
8ICR3	Troponin IC complex, cardiac, human	Recombinant	>95%
8IFC20	Troponin I (fragment 28-110) – troponin C complex, cardiac, human	Recombinant	>95%
8ITCR	Troponin ITC complex, cardiac, human	Recombinant	>95%

# Soluble CD40 ligand (sCD40L)

CD40 ligand (CD40L) is a member of the tumor necrosis factor (TNF) family and is expressed on the surface of CD4+ T-cells, basophiles, platelets and mast cells. The binding of CD40L to its receptor CD40 mediates various inflammatory processes in cells. Soluble CD40 ligand (sCD40L) is formed upon the cleavage of CD40L. This soluble form has been shown to act as a cytokine. Some studies suggest that an increased level of sCD40L in blood might correlate with acute and chronic heart failure and the severity of the disease.

## MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4CD40*	Monoclonal mouse anti-soluble CD40 ligand (sCD40L)	Enzyme immunoassays Western blotting

\* Several MAbs available under one catalogue number. Please see [www.hytest.fi](http://www.hytest.fi).

# Glycogen phosphorylase isoenzyme BB (GPBB)

Glycogen phosphorylase isoenzyme BB (GPBB) plays an important role in glycogen turnover. BB isoform is synthesized by cardiac and brain tissues. GPBB is considered to be an early marker of myocardial cell death and its release kinetics closely resemble those of myoglobin and FABP. GPBB could be a useful marker in the diagnosis of acute coronary syndrome and unstable angina.

## MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4GP31	Monoclonal mouse anti-glycogen phosphorylase isoenzyme BB (GPBB)	Enzyme immunoassays Western blotting

# Soluble lectin-like oxidized LDL receptor (sLOX-1)

Soluble lectin-like oxidized LDL receptor (sLOX-1) is produced by the proteolytic cleavage of the extracellular domain of LOX-1. LOX-1 is a transmembrane protein that is for example found on the cell surface of endothelial cells and smooth muscle cells. The serum level of sLOX-1 is increased in atherosclerotic conditions and with inflammation. A few studies suggest that sLOX-1 could serve as a biomarker for plaque rupture and that it might have clinical value in diagnosing atherosclerosis-related diseases.

## MONOCLONAL ANTIBODIES

Cat.#	Product	Tested applications
4LOX1	Monoclonal mouse anti-sLOX-1	Enzyme immunoassays Western blotting

## REFERENCES

1. **Alhadi HA and Fox KA.** Do we need additional markers of myocyte necrosis: the potential value of heart fatty-acid-binding protein, *QJMed.* 2004; 97: 187–198.
2. **Ballantyne, C.M., et al.** Lipoprotein-associated phospholipase A2, high-sensitivity C-reactive protein, and risk for incident coronary heart disease in middle-aged men and women in the Atherosclerosis Risk in Communities (ARIC) study. *Circulation,* 2004, 109(7):837-842.
3. **Bayes-Genis A. et al.** Pregnancy-associated plasma protein A as a marker of acute coronary syndromes.// *N Engl J Med.* 2001, 345(14): 1022–9.
4. **Body R and Ferguson C.** Pregnancy-associated plasma protein A: a novel cardiac marker with promise. *Emerg Med J.* 2006, 23(11): 875–877.
5. **Brennan ML et al.** Prognostic value of myeloperoxidase in patients with chest pain. *New Eng J Med.* 2003, 349(17): 1595–1604.
6. **Brilakis, E.S., et al.** Association of lipoprotein-associated phospholipase A2 levels with coronary artery disease risk factors, angiographic coronary artery disease, and major adverse events at follow-up. *Eur Heart J,* 2005, 26(2):137-144.
7. **Clearfield MB.** C-reactive protein: a new risk assessment tool for cardiovascular disease. *JAOA* 2005; 105(9):409–416.
8. **Garg, P.K., et al.** Lipoprotein-Associated Phospholipase A2 and Incident Peripheral Arterial Disease in Older Adults: The Cardiovascular Health Study. *Arterioscler Thromb Vasc Biol,* 2016, 36(4):750-756.
9. **Hájek P et al.** Pregnancy-associated plasma protein A and proform eosinophilic major basic protein in the detection of different types of coronary artery disease. *Physiol. Res.* 2008, 57(1): 23–32.
10. **Heeschen C et al.** Pregnancy-associated plasma protein-A levels in patients with acute coronary syndromes Comparison with markers of systemic inflammation, platelet activation, and myocardial necrosis. *J. Am. Coll. Cardiol.* 2005, 45:229–237.
11. **Katrukha AG.** Antibody selection strategies in cardiac troponin assays. *Cardiac Markers,* 2nd edition, Edited by Alan HB. Wu. 2003, 173–185.
12. **Katrukha AG et al.** Troponin I is released in bloodstream of patients with acute myocardial infarction not in free form but as complex. *Clin. Chem.* 1997, 43(8): 1379-1385.
13. **Katrukha AG et al.** Degradation of cardiac troponin I: implication for reliable immunodetection. *Clin. Chem.* 1998, 44(12): 2433-2440.
14. **Koenig W et al.** Prospective Study of High-Sensitivity C-Reactive Protein as a Determinant of Mortality: Results from the MONICA/KORA Augsburg Cohort Study, 1984–1998. *Clin. Chem.* 2008; 54:2:335–342.
15. **Meyer MH et al.** CRP determination based on a novel magnetic biosensor. *Biosens. Bioelectron.* 2007, 22(6):973–979.
16. **Nambi V.** The use of myeloperoxidase as a risk marker for atherosclerosis. *Curr. Atheroscler. Rep.* 2005, 7(2):127–131.
17. **Penttilä K et al.** Myoglobin, creatine kinase MB isoforms and creatine kinase MB mass in early diagnosis of myocardial infarction in patients with acute chest pain. *Clin. Biochem.* 2002, 35(8):647–653.
18. **Postnikov AB et al.** N-terminal and C-terminal fragments of IGFBP-4 as novel biomarkers for short-term risk assessment of major adverse cardiac events in patients presenting with ischemia. *Clin. Chem.* 2012, 45:519–524.
19. **Qin QP et al.** Release Patterns of Pregnancy Associated Plasma Protein A (PAPP-A) in Patients with Acute Coronary Syndromes. *Scand Cardiovasc. J.* 2002, 36(6):358–361.
20. **Ridker PM.** C-reactive protein: a simple test to help predict risk of heart attack and stroke. *Circulation.* 2003. 108:e81–e85.
21. **Scirica BM et al.** Clinical Application of C-Reactive Protein Across the Spectrum of Acute Coronary Syndromes. *Clin. Chem.* 2007, 53:1800–1807.
22. **Seferian KR et al.** Immunodetection of glycosylated NT-proBNP circulating in human blood. *Clin. Chem.* 2008, 54(5):866–873.
23. **Shiesh SC et al.** Determination of C-reactive protein with an ultra-sensitivity immunochemiluminometric assay. *J. Immunol. Methods.* 2006, 311(1–2):87–95.
24. **Sin KK et al.** Fluorogenic nanocrystals for highly sensitive detection of C-reactive protein. *IEE Proc. Nanobiotechnol.* 2006, 153(3):54–58.
25. **Stafforini, D.M.** Biology of platelet-activating factor acetylhydrolase (PAF-AH, lipoprotein associated phospholipase A2). *Cardiovasc Drugs Ther,* 2009, 23(1):73-83.
26. **Stafforini, D.M.** Plasma PAF-AH (PLA2G7): Biochemical Properties, Association with LDLs and HDLs, and Regulation of Expression. *Enzymes,* 2015, 38:71-93.
27. **Storch J and Thumser AE.** The fatty acid transport function of fatty acid-binding proteins. *Biochim. Biophys. Acta* 2000; 1486: 28–44
28. **Terkelsen CJ. et al.** Temporal course of pregnancy-associated plasma protein-A in angioplasty-treated ST-elevation myocardial infarction patients and potential significance of concomitant heparin administration. *Am. J. Cardiol.* 2009, 103(1):29–35.
29. **Tertti R et al.** Intravenous administration of low molecular weight and unfractionated heparin elicits a rapid increase in serum pregnancy-associated plasma protein A. *Clin Chem.* 2009, 55(6):1214–1217.
30. **van Nieuwenhoven FA et al.** Discrimination between myocardial and skeletal muscle injury by assessment of the plasma ratio of myoglobin over fatty acid-binding protein. *Circulation.* 1995, 92(10):2848–2854.
31. **Zhang R. et al.** Association between myeloperoxidase levels and risk of coronary artery disease. *JAMA.* 2001, 286(17):2136–2142.

# Selected cardiac marker articles from Hytest scientists

**Filatov VL, et al.** Epitope mapping of anti-troponin I monoclonal antibodies. *Biochem. Mol. Biol. Int.* 1998, 45(6): 1179-1187.

**Filatov VL, et al.** Troponin: structure, properties, and mechanism of functioning. *Biochemistry.* 1999, 64(9): 969-985.

**Katrakha A, et al.** Biochemical factors influencing measurement of cardiac troponin I in serum. *Clin. Chem. Lab. Med.* 1999, 37(11-12): 1091-1095.

**Katrakha A, et al.** New approach to standardization of human cardiac troponin I (cTnI). *Scand. J. Clin. Lab. Invest., Suppl.* 1999, 230: 124-127.

**Katrakha AG, et al.** A new method of human cardiac troponin I and troponin T purification. *Biochem. Mol. Biol. Int.* 1995, 36: 195-202.

**Katrakha AG, et al.** Troponin I is released in bloodstream of patients with acute myocardial infarction not in free form but as complex. *Clin. Chem.* 1997, 43(8): 1379-1385.

**Katrakha AG, et al.** Degradation of cardiac troponin I: implication for reliable immunodetection. *Clin. Chem.* 1998, 44(12): 2433-2440.

**Katrakha AG.** Antibody selection strategies in cardiac troponin assays. *Cardiac Markers*, 2nd edition, Edited by Alan H. B. Wu. 2003, 173-185.

**Katrakha IA, et al.** Thrombin-Mediated Degradation of Human Cardiac Troponin T. *Clin Chem.* 2017, 63(6):1094-1100.

**Konev AA, et al.** Characterization of endogenously circulating IGFBP-4 fragments—Novel biomarkers for cardiac risk assessment. *Clin. Biochem.* 2015, 48(12): 774-780.

**Postnikov AB, et al.** N-terminal and C-terminal fragments of IGFBP-4 as novel biomarkers for short-term risk assessment of major adverse cardiac events in patients presenting with ischemia. *Clin. Chem.* 2012, 45: 519-524.

**Røsjo H, et al.** Diagnostic utility of a single-epitope sandwich B-type natriuretic peptide assay in stable coronary artery disease: data from the Akershus Cardiac Examination (ACE) 1 Study. *Clin. Biochem.* 2012 45(16-17): 1269-1275.

**Schulz O, et al.** Clinical differences between total PAPP-A and measurements specific for the products of free PAPP-A activity in patients with stable cardiovascular disease. *Clin. Biochem.* 2013, 46; doi: 10.1016/j.clinbiochem.2013.10.027

**Seferian KR, et al.** The brain natriuretic peptide (BNP) precursor is the major immunoreactive form of BNP in patients with heart failure. *Clin. Chem.* 2007, 53: 866-873.

**Seferian KR, et al.** Immunodetection of glycosylated NT-proBNP circulating in human blood. *Clin. Chem.* 2008, 54(5): 866-873.

**Semenov AG and Katrukha AG.** Different Susceptibility of B-Type Natriuretic Peptide (BNP) and BNP Precursor (proBNP) to Cleavage by Neprilysin: The N-Terminal Part Does Matter. *Clin Chem.* 2016, 62(4):617-622.

**Semenov AG and Seferian KR.** Biochemistry of the human B-type natriuretic peptide precursor and molecular aspects of its processing. *Clin. Chim. Acta.* 2011, 412(11-12): 850-860.

**Semenov AG, et al.** Processing of pro-brain natriuretic peptide is suppressed by O-glycosylation in the region close to the cleavage site. *Clin. Chem.* 2009, 55(3): 489-498.

**Semenov AG, et al.** Processing of Pro-B-Type Natriuretic Peptide: Furin and Corin as Candidate Convertases. *Clin. Chem.* 2010, 56(7): 1166-1176.

**Semenov AG, et al.** Human Pro-B-Type Natriuretic Peptide Is Processed in the Circulation in a Rat Model. *Clin. Chem.* 2011, 57(6): 883-890.

**Semenov AG, et al.** Searching for a BNP standard: Glycosylated proBNP as a common calibrator enables improved comparability of commercial BNP immunoassays. *Clin Biochem.* 2017, 50(4-5):181-185.

**Tamm NN, et al.** (2008) Novel immunoassay for quantification of brain natriuretic peptide and its precursor in human blood. *Clin. Chem.* 54(9), 1511-1518.

**Tamm NN, et al.** Measurement of B-type natriuretic peptide by two assays utilizing antibodies with different epitope specificity. *Clin. Biochem.* 2011, 44(2-3): 257-259.

**Vylegzhanina AV, et al.** Epitope Specificity of Anti-Cardiac Troponin I Monoclonal Antibody 8I-7. *Clin. Chem.* 2013, 59(12): 1814-1816.

**Vylegzhanina AV, et al.** Anti-Cardiac Troponin Autoantibodies Are Specific to the Conformational Epitopes Formed by Cardiac Troponin I and Troponin T in the Ternary Troponin Complex. *Clin. Chem.* 2017, 63(1), 343-350.

# Patents and trademarks

Immunoassay Kit for Quantification of BNP and proBNP (US 9,145,459)

Detection of Cardiac Muscle Necrosis by Immunoassay and Appropriate Antibodies (EP 0965043, FI 104857)


Method and Kit for the Diagnosis of Troponin I (US 7,285,418, EP 0938678)

Stable Standards for BNP Immunoassays (EP 2084544, CN 101641601, CA 2669024)

Immunoassay for Quantification of an Unstable Antigen Selected from BNP and proBNP (US 9,034,591, US 9,034,592, JP 5686593, CN 101842707, CA 268391, EP 2135087)

Detection of IGFBP-4 Fragments as a Diagnostic Method (EP 2448969, US 9,012,610, JP 5840605)

SES-BNP™ is a trademark or registered trademark of Hytest Ltd.



# More information

More detailed information regarding the performance of our products, a full list of individual MAbs and recommendations for capture-detection antibody pairs (when available) can be found on our website — [www.hytest.fi](http://www.hytest.fi).

**You are also most welcome to contact our Sales Team directly by writing to [hytest@hytest.fi](mailto:hytest@hytest.fi).**





# Scientific excellence for IVD

#### **HYTEST LTD**

Intelligate 1, 6th floor, Joukahaisenkatu 6  
FI-20520 Turku, FINLAND  
Tel. +358-2-512-0900 • E-mail: [hytest@hytest.fi](mailto:hytest@hytest.fi)  
**HYTEST.FI**

#### **HYTEST NORTH AMERICA INC.**

800 MacArthur Blvd, STE 100,  
Mahwah, New Jersey 07430  
Tel. 833-368-3489 • E-mail: [hytestna@hytest.fi](mailto:hytestna@hytest.fi)  
**HYTESTNA.COM**

#### **OUR BROCHURES**

Product Catalog 2026 • Cardiac Markers •  
Neuroscience • Inflammation •  
Infectious Diseases • Metabolic  
Syndrome • Thyroid Diseases • Veterinary