



# Welcome to HyTest Webinar COVID-19: Biomarkers to distinguish severity and prognosis

May 7, 2020

C O N F I D E N T I A L


- 1 Introductions: Company and Speaker *by Dr. Netta Fatal*
- 2 Talk 45 min *by Dr. Alexander Semenov*
- 3 Q&A 15 min *by Dr. Alexander Semenov and Dr. Alexander Postnikov*

# HyTest Ltd.

*One of the key  
raw material  
suppliers for the  
IVD industry*

We develop and produce monoclonal antibodies and antigens that are mainly used as key components in laboratory tests.

HyTest was established in 1994.

-  Headquarters in Finland, operations in China, Russia and North America
-  Sales to over 50 countries
-  Active participation in IFCC and AACC standardization committee work
-  Operations ISO 9001:2015 compliant

# Comprehensive product line

## PRODUCT CATEGORIES

- Monoclonal antibodies
- Polyclonal antibodies
- Antigens
- Plasma and serum
- Over 1,000 different reagents

## KEY PRODUCT AREAS



Cardiac  
Markers



Metabolic  
Syndrome



Infectious  
Diseases



Inflammation



Veterinary

## OTHER PRODUCT AREAS



Blood Coagulation  
and Anemia



Immunology  
and Serology



Fertility and  
Pregnancy



Hormones



Tumor Markers



Neuroscience



Gangliosides



Kidney diseases



Inflammation



Microbial and  
Plant Toxins



Biodefence



Molecular Biology

# Today's speaker: Dr. Alexander Semenov



- MSc. and PhD degrees from the Moscow State University (MSU)
- Senior Scientist and Project Manager, joined HyTest R&D in 2005
- Involved in research projects focused on BNP and NT-proBNP, in charge of developing antibodies and immunoassays specific to these heart failure biomarkers
- Co-author of about 20 publications and patents
- Member of The Joint Committee for Traceability in Laboratory Medicine (JCTLM) Proteins Review Team

# Panelist Dr. Alexander Postnikov



- MSc and PhD degrees from the Moscow State University (MSU)
- Joined HyTest R&D in 2006, manages the group of development and production of recombinant antibodies
- Involved in several NPD projects including invention of novel biomarkers for cardiac risk assessment (IGFBP-4 fragments)
- Co-author of about 30 scientific publications and patents
- Member of the IFCC Working Group on PAPP-A Standardization



# COVID-19: Biomarkers to distinguish severity and prognosis


Alexander G. Semenov, PhD

CONFIDENTIAL



# Coronavirus disease 2019 (COVID-19) pandemic

- In December 2019, a cluster of atypical pneumonia patients epidemiologically linked to a wholesale market in Wuhan (Hubei Province, China) was detected
- Initially, a novel coronavirus was called 2019-nCoV
- Later it was termed the SARS-CoV-2 virus as it is very similar to the one that caused the outbreak of severe respiratory disease (SARS) in 2003





# Coronavirus disease 2019 (COVID-19) pandemic

- At the end of January the World Health Organization (WHO) declared the new infectious disease COVID-19 a global health emergency
- On 11 March 2020 the new infectious disease was recognized as a pandemic by the WHO
- In up to 15% of infected patients the clinical course may be complicated by the onset of a severe form of interstitial pneumonia, which may then progress towards acute respiratory distress syndrome and/or multi organ failure and death



1

SARS-CoV-2 virus

2

Nucleic acid testing and serological assays for SARS-CoV-2

3

Effects of SARS-CoV-2 on the cardiovascular and other systems

4

COVID-19 and cardiac biomarkers

5

COVID-19 and inflammatory biomarkers

6

COVID-19 and coagulation, kidney and muscle injury biomarkers

1

## SARS-CoV-2 virus

2

Nucleic acid testing and serological assays for SARS-CoV-2

3

Effects of SARS-CoV-2 on the cardiovascular and other systems

4

COVID-19 and cardiac biomarkers

5


COVID-19 and inflammatory biomarkers

6

COVID-19 and coagulation, kidney and muscle injury biomarkers


# SARS-CoV-2 is a novel coronavirus of zoonotic origin

- SARS-CoV-2 belongs to a large family of single-stranded RNA viruses (+ssRNA)
- It has a crown-like appearance under an electron microscope (*coronam* is the Latin term for crown) due to the presence of spike glycoproteins on the envelope
- Beta-coronaviruses (like SARS-CoVs) can cross species barriers and can cause in humans illness ranging from the common cold to more severe diseases such as SARS (2003) and MERS (2012)



# Mechanisms of entry of SARS-CoV-2

- SARS-CoV-2 enters the body through the nose and throat (e.g. when virus-laden droplets are inhaled)
- SARS-CoV-2 enters cells through an interaction with angiotensin-converting enzyme 2 (ACE-2)
- ACE-2 is a transmembrane peptidase that degrades angiotensin II to generate angiotensin 1-7 (vasodilator)



# Should inhibitors of the renin–angiotensin system be withdrawn in patients with COVID-19?

- ACE and ACE-2 belong to the same peptidase family, however, they have two very different physiological functions. ACE is a target for the treatment of hypertension.
- It's hypothesized that ACE inhibitors (ACE-Is) could act as a potential risk factor for COVID-19 by up-regulating ACE-2
- Patients taking ACE-Is or ARBs (angiotensin-receptor blockers) may be more susceptible for viral infection and have higher mortality because they are older, more frequently hypertensive, diabetic, and/or having renal disease
- It's suggested that ACE-I and ARB therapy should be maintained or initiated in patients with heart failure, hypertension, or myocardial infarction, according to current guidelines, irrespective of SARS-CoV-2



1

SARS-CoV-2 virus

2

Nucleic acid testing and serological assays for SARS-CoV-2

3

Effects of SARS-CoV-2 on the cardiovascular and other systems

4

COVID-19 and cardiac biomarkers

5

COVID-19 and inflammatory biomarkers

6


COVID-19 and coagulation, kidney and muscle injury biomarkers

# Assays available for the detection of SARS-CoV-2

- Nucleic acid testing:
  - The most sensitive method combined with high specificity and high efficiency
  - The limit of detection reaches 100 copies/mL
  - Stringent performance assessment is still an urgent need
- Serological testing:
  - Testing for specific IgM, IgG or viral antigens (ELISA, CLIA, etc. and rapid serological testing)
  - Help clinically discriminate among infections when the nucleic acid testing result is negative
  - Sensitive and specific serological assays are not as easily established as nucleic acid testing assays

# Nucleic acid testing for the detection of SARS-CoV-2

- The assays used in many laboratories are real time PCR assays targeting two different amplification regions:
  - The E (envelope protein)
  - RdRp (RNA-dependent RNA polymerase) genes
- The FDA approved assay from Abbott Diagnostics uses an isothermal nucleic acid amplification method



# Nucleic acid testing for the detection of SARS-CoV-2

- Viral load peaks in the first week of disease onset
- Viral RNA can be detected in patients in the 2nd week of disease onset, but the viral load is low
- The optimal specimen type for SARS-CoV-2 detection is yet to be determined
- Saliva samples may be better than nasal/throat swabs, however, additional studies are needed

# Limitations of nucleic acid testing for Covid-19


- Require high-quality viral RNA in sufficient amounts (amounts vary tremendously between patients)
- Viral RNA is sensitive to spoilage during collection, transport, and storage
- Despite high sensitivity, a negative nucleic acid testing is insufficient to exclude SARS-CoV-2 infection in patients with high clinical suspicion
- If a negative nucleic acid testing is observed at one or two time points, other approaches for testing should be considered, including specific IgM and IgG assays (serological assays)

# Serological testing of SARS-CoV-2


- Serology tests are used for the detection of antibodies in the blood in order to diagnose an active or previous infection
- These antibodies can be detected only several days after the exposure to the disease causing agent (e.g. SARS-CoV-2)
- Tests are suitable for assessing whether the person has been exposed to the virus at some point and developed immune response against it



# Principle of serological testing of SARS-CoV-2



Detection of anti-SARS-CoV-2 antibodies in serum samples enables to confirm the immune response.




# Serological testing of SARS-CoV-2

- In general, IgM antibodies can be detected earlier than IgG antibodies, however, studies related to SARS-Cov-2 show mixing results
- Serology tests are useful for:
  - Determining how widely a disease has spread within a community
  - Determining whether the person already has developed immune response to the disease
  - Determining who may donate their blood, which may serve as a possible treatment for those who are seriously ill from COVID-19
  - Monitoring the effectiveness of vaccine (ability to raise immune response)

# Immunoassays to detect SARS-CoV-2 antigen

- Rapid antigen tests may provide the advantage of fast time to results and low-cost detection, but require a high viral count
- Immunoassays for SARS-CoV-2 antigen are still new to the market
- The viral spike protein can be used as an antigen
- Nucleocapsid protein is a promising target for clinical diagnostics:
  - Mutations are much less common than in spike protein
  - Produced at high levels within infected cells

# Levels of SARS-CoV-2 RNA/antigen, IgM and IgG after infection



1

SARS-CoV-2 virus

2

Nucleic acid testing and serological assays for SARS-CoV-2

3

Effects of SARS-CoV-2 on the cardiovascular and other systems

4

COVID-19 and cardiac biomarkers

5

COVID-19 and inflammatory biomarkers

6

COVID-19 and coagulation, kidney and muscle injury biomarkers

# Effects of SARS-CoV-2 on different organs

*Strokes, seizures,  
confusion, and brain  
inflammation*

*Inflamed and  
damaged alveolus*

*Infection of  
gastrointestinal  
tract causes  
diarrhea*



*Conjunctivitis/  
inflammation*

*Blood clots, heart  
attacks, and cardiac  
inflammation*

*Kidney damage*

*Liver disfunction*



# Effects of SARS-CoV-2 on the cardiovascular system

- Acute respiratory infections are well-recognized triggers for cardiovascular diseases (CVD), and the underlying CVD is usually associated with comorbidities, which may increase the incidence and severity of infectious diseases.
- Advanced age (>60 years), male sex, and presence of comorbidities are known to be the major risk factors for COVID-19 mortality
- Presence of cardiac injury (defined by elevated troponin levels), myocarditis, and ARDS are other strong and independent factors associated with mortality.

# Effects of SARS-CoV-2 on the cardiovascular system

- Heart damage and arrhythmias are common in patients hospitalized for COVID-19
- Abnormal blood clotting (causing pulmonary embolism or stroke) and blood vessel constriction have been reported
- How SARS-CoV-2 attacks the heart and blood vessels is unknown:
  - Heart and blood vessels are rich in ACE-2 (the target of SARS-CoV-2)
  - Lack of oxygen may damage blood vessels
  - Cytokine storm may cause damage to the heart

1

SARS-CoV-2 virus

2

Nucleic acid testing and serological assays for SARS-CoV-2

3

Effects of SARS-CoV-2 on the cardiovascular and other systems

4

COVID-19 and cardiac biomarkers

5

COVID-19 and inflammatory biomarkers

6

COVID-19 and coagulation, kidney and muscle injury biomarkers


# Myocardial injury is a common comorbidity in patients with COVID-19

- Rise and/or fall of troponin indicating myocardial injury is common among patients with acute respiratory infections and correlated with disease severity
- Abnormal troponin values are common among those with COVID-19 infection particularly when testing with a high sensitivity cardiac troponin (hs-cTn) assay
- Detectable hs-cTnI was observed in most patients with COVID-19, and hs-cTnI was significantly elevated in more than half of the patients that died

*Januzzi J. ACC, 2020.*

# Cardiac troponin I (cTnI) levels in patients with COVID-19


- cTnI values are significantly increased in patients with severe SARS-CoV-2 infection compared to those with milder forms of disease
- Adjunctive cardioprotective therapies may be advisable in patients with significant elevation of cardiac injury biomarkers



Meta-analysis by Lippi G. et al. *Progress in Cardiovascular Diseases*, 2020.

# Cardiac troponin T (cTnT) testing in patients with COVID-19

- Myocardial injury was significantly associated with fatal outcome of COVID-19, while the prognosis of patients with underlying CVD, but without myocardial injury was relatively favorable
- Aggressive treatment may be considered for patients at high risk of myocardial injury



Guo T. et al. JAMA Cardiology, 2020.



- *Data from Wuhan, China, confirm that cardiac injury is common in patients hospitalized with COVID-19 and, unsurprisingly, is tightly linked to higher mortality.*
- *Initial measurement of cardiac damage biomarkers immediately after hospitalization for SARS-CoV-2 infection may help identifying a subset of patients with possible cardiac injury.*
- *Monitoring of cTnI/cTnT during hospital stay may also be helpful to predict the progression of COVID-19 towards a worse clinical picture.*

*High-sensitivity cardiac troponin “could actually be an ally in the fight against COVID-19”.*


*Troponin is a “crucial diagnostic and prognostic aid in what will become even more challenging times for healthcare provision worldwide”.*

**Andrew R Chapman**

Cardiologist in BHF Centre for Cardiovascular  
Science, Royal Infirmary of Edinburgh, Edinburgh, UK

# Natriuretic peptide response in patients with COVID-19


- Natriuretic peptides are biomarkers of myocardial stress and are frequently elevated among patients with severe respiratory illnesses typically in the absence of elevated filling pressures or clinical heart failure
- Much like troponin, elevation of BNP or NT-proBNP is associated with an unfavorable course among patients with ARDS



# NT-proBNP levels in patients with COVID-19


- Patients with higher NT-proBNP (above 88.64 pg/mL) level had more risks of in-hospital death
- After adjusting for potential cofounders in separate modes, NT-proBNP presented as an independent risk factor of in-hospital death in patients with severe COVID-19

*Gao L. et al. Medrxiv, 2020.*



# NT-proBNP levels in patients with COVID-19

- The cutoff value of NT-proBNP to predict the adverse outcome of severe COVID-19 patients was far lower than the threshold to diagnose heart failure
- ***The prognostic effect of NT-proBNP might indicate the extent of cardiac stress and inflammation***
- NT-proBNP level might be helpful to identifying patients with poor prognoses early



Gao L. et al. Medrxiv, 2020.

# BNP levels in patients with COVID-19

- Patients with COVID-19 often have abnormal BNP/NT-proBNP in plasma
- Compared with patients with normal BNP, patients with high BNP (>100 pg/mL) were more likely to develop severe pneumonia, invasive mechanical ventilation, continuous renal replacement therapy, and be admitted to the intensive care unit
- ***The level of BNP in plasma may reflect the severity of inflammation and stress.*** This may partly explain why patients with high plasma BNP levels had a bad outcomes

*Liu Y. et al. Medrxiv, 2020.*

1

SARS-CoV-2 virus

2

Nucleic acid testing and serological assays for SARS-CoV-2

3

Effects of SARS-CoV-2 on the cardiovascular and other systems

4

COVID-19 and cardiac biomarkers

5

COVID-19 and inflammatory biomarkers

6

COVID-19 and coagulation, kidney and muscle injury biomarkers

# Inflammatory biomarkers

- C-reactive protein (CRP) is routinely used as a non-specific marker of inflammation
  - To be distinguished from high-sensitivity CRP (hsCRP) that used as a marker of increased risk for cardiac diseases
- Procalcitonin (PCT) is a marker of disorders that are accompanied by systemic inflammation and sepsis
- Serum amyloid A (SAA) proteins form a family of apolipoproteins. Some SAAs are expressed constitutively, while others are expressed in response to inflammation




# C-reactive protein (CRP) testing in patients with COVID-19

- CRP does not normally elevate significantly in mild viral respiratory infections
- However, significant increase of CRP has also been reported in COVID-19 patients
- CRP testing may be useful in the initial evaluation of coronavirus patients
- CRP testing and complete blood count can be also used at the forefront for evaluation of infection and to direct patients further on the treatment path


*Zhang et al. Lancet Respir Med, 2020.*

Pentameric  
plasma CRP



# Procalcitonin (PCT) testing in patients with COVID-19

- PCT values are not substantially modified in patients with viral infections
- Serial PCT measurement in patients with or without severe COVID-19 may play a role for predicting evolution towards a more severe form of disease
- Evidence from a meta-analysis



*Lippi G. et al. Clin Chem Acta, 2020.*

# CRP, PCT and ferritin testing in patients with COVID-19

- The serum levels of CRP, PCT and ferritin are markedly increased in very severe compared with severe COVID-19
- Increased CRP, PCT and ferritin level might correlate to secondary bacterial infection and associated with poor clinical prognosis

*Zhou B. et al. Research Square, 2020.*

# SAA is a biomarker to distinguish the severity and prognosis of COVID-19

- Patients with higher initial SAA are more likely to have poor CT imaging
- The level of SAA and CRP significantly increased in patients with COVID-19
- As disease progressed from mild to critically severe, SAA and CRP gradually increased
- SAA and lymphocytes are valuable in predicting the severity and distinguishing critically ill patients from mild ones

*Li H. et al Journal of infection, 2020.*

1

SARS-CoV-2 virus

2

Nucleic acid testing and serological assays for SARS-CoV-2

3

Effects of SARS-CoV-2 on the cardiovascular and other systems

4

COVID-19 and cardiac biomarkers

5


COVID-19 and inflammatory biomarkers

6

COVID-19 and coagulation, kidney, and muscle injury biomarkers

# D-dimer levels in patients with COVID-19

- D-dimer is a biomarker of pathological coagulation that underlies the pathogenesis of most cardiovascular diseases
- D-dimer values are frequently enhanced in patients with COVID-19, being variably observed in 36 to 43% of positive cases
- D-dimer values are considerably higher in COVID-19 patients with severe disease than in those without the disease



*Lippi G. et al. Thrombosis and Haemostasis, 2020.*

# D-dimer levels in patients with COVID-19

- D-dimer measurement may be helpful in predicting evolution toward worse clinical picture in COVID-19 patients
- D-dimer measurement may help to define whether adjunctive antithrombotic therapies (e.g., anticoagulants, antithrombin or thrombomodulin) might be helpful in patients with severe COVID-19

*Lippi G. et al. Thrombosis and Haemostasis, 2020.*

---

*“The more we look, the more likely it becomes that blood clots are a major player in the disease severity and mortality from COVID-19.”*

**Behnood Bikdeli**

Cardiovascular medicine fellow at Columbia University  
Medical Center

*Meredith W. et al. Science, 2020*



# Cystatin C levels in patients with COVID-19

- Patients infected with COVID-19 has renal failure is a frequent complication of acute respiratory distress syndrome
- Serum biomarkers, including urea, creatinin, cystatin C, which reflect glomerular filtration function, may be useful as potential indicators for the early diagnosis of severe COVID-19 and to distinguish it from mild COVID-19

*Xiang J. et al. MedRxiv, 2020.*

# Elevated levels of myoglobin in association with COVID-19 severity


- Several studies have reported elevated levels of myoglobin in association with COVID-19 severity
- SARS-CoV-2 infection might induce a myositis similar to that observed in severe influenza infections
- Rapid clinical recognition of muscle injury in patients with COVID-19 can be lifesaving

*Bangash M. et al. Lancet Gastroenterol Hepatol, 2020.*

*Jin M. et al. Emerging Infectious Diseases, 2020.*


# Albumin levels in patients with COVID-19

- Serum albumin (HSA) levels were associated with an increased risk of death
- Decreased level of HSA may be a medical sign of decreased production in the liver, increased loss in the gastrointestinal tract or kidneys, increased use in the body
- The combinations of the hypo-albuminemia, lymphopenia, and high concentrations of CRP and LDH in SARS-CoV-2 infected patients upon hospital admission may predict more severe acute lung injury



Liu Y. et al. *Sci China Life Sci*, 2020.

# COVID-19: biomarkers to distinguish severity and prognosis



# Coronavirus as if painted by famous artists



Van Gogh



Picasso



Dali



de Vinci



*Ceci n'est pas un Coronavirus*


Magritte



Pollock



Warhol



Mondrian

VALOTT

*Similarly, the response to the virus is totally different in different individuals.*


*Biomarkers can be used as precise tools to distinguish individual differences and guide therapy accordingly.*

# Concluding remarks

- Patients with comorbidities (like diabetes mellitus, hypertension, cardiovascular, chronic lung and chronic kidney disease) are particularly susceptible to COVID-19 infection and are likely to have more severe illness
- COVID-19 is associated with multiple direct and indirect cardiovascular complications including acute myocardial injury, myocarditis, arrhythmias and venous thromboembolism
- Future therapies for COVID-19 may have cardiovascular and other systemic side effects

# Concluding remarks

- The key role of biomarkers may be to predict the progression of COVID-19 towards a worse clinical picture and provide a signal to manage patients differently
- They can help to triage patients to critical care, guide the use of supportive treatments, and facilitate targeted investigations in those most likely to benefit
- Careful monitoring of biomarkers is important in reducing the complications and mortality in COVID-19 patients



# Products available from HyTest

- For serology assays:
  - Antibodies specific to IgM and IgG
  - SARS-CoV-2 antigens **in development** (ask for availability)
- SARS-Cov-2 specific antibodies **in development**
- Antibodies and native/recombinant antigens for immunoassay development covering several clinical categories:
  - Cardiac markers
  - Inflammatory markers
  - Blood coagulation
  - Infectious diseases
  - Etc.





From the science  
in our hearts to  
a healthier world.  
Together with you.  
Today and tomorrow.

[www.hytest.fi](http://www.hytest.fi)  
[hytest@hytest.fi](mailto:hytest@hytest.fi)

