



DataSheet

CATALOGUE #: PSN5

PRODUCT NAME: Polyclonal anti-SARS-CoV-2 Nucleoprotein – TEST SAMPLE

Immunogen: Recombinant SARS-CoV-2 Nucleoprotein

Host animal: Goat

Specificity: SARS-CoV-2 Nucleoprotein

Applications: Immunoassays. Recommended pairs for sandwich immunoassay (capture-detection):

PS5N – C518

PS5N – C527

(C518 and C527 are under Cat.# 3CV4)

Purification: Immunoaffinity method utilizing SARS-CoV-2 Nucleoprotein-sepharose as an affinity matrix

Presentation: PBS, pH 7.4, 0.09 % sodium azide (NaN₃)

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

www.hytest.fi | hytest@hytest.fi