



DataSheet

CATALOGUE #: 8COV5

PRODUCT NAME: SARS-CoV-2 nucleoprotein fragment N47-A173, recombinant – TEST SAMPLE

Source: Expressed in *E. coli*. Contains His₆-tag.

Presentation: 8 M urea, 250 mM NaCl, and 250 mM of imidazole, pH 10

Storage: -70°C (-65 ... -80 °C allowed)

Other information: Avoid repeated freezing and thawing. It is recommended to aliquot the product after the first thawing.

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.